



LFH-W Series (Walk-in)


Available: 39", 50", 62"



- With memory function in case of power-failure.
- Operators are available to walk into the fume hood to operate if need.
- Large room exhaust system, suitable for large laboratory instruments and equipment.



- Water Sink



- Water & Gas Remote Control Valve



- Application



- Water & Gas Tap



Walk-in Fume Hood

Model		LFH-4W	LFH-5W	LFH-6W
External Size (WxDxH)		1200×800×2250mm 47.2"×31.5"×88.6"	1500×800×2250mm 59"×31.5"×88.6"	1800×800×2250mm 70.9"×31.5"×88.6"
Internal Size (WxDxH)		980×550×1670mm 38.6"×21.6"×65.7"	1280×550×1670mm 50.4"×21.6"×65.7"	1580×550×1670mm 62.2"×21.6"×65.7"
Max Opening		1300mm (51.2")		
Exhaust Airflow Volume		1150m ³ /h	1500m ³ /h	1850m ³ /h
Inflow Air Velocity		0.3~0.8m/s		
Noise		≤62dB		
Vibration		Amplitude≤5μm		
Illumination		≥800Lux		
LED Lamp		12W×2	16W×2	16W×2
Power Supply		AC220V±10%, 50/60HZ; 110V±10%, 60Hz		
Consumption		400W	500W	800W
Front Window		Manual, 5mm toughened glass, height adjustable		
Material	Exterior	1.0mm cold-rolled steel with bacteria powder coating		
	Interior	Compact laminate board with good acid and alkali resistance function		
Standard Accessory		External blower, LED lamp, Water tap, Gas tap, PP cup sink, Remote control water, valve, Remote control gas valve, Waterproof sockets*4, 4.3 meters Ø300mm PVC exhaust duct		
Sockets		Total load of 4 waterproof sockets; ≤500W		
Package Size (WxDxH)	Main Body	1400×1050×2450mm 55.1"×41.3"×96.5"	1710×1050×2450mm 67.3"×41.3"×96.5"	2010×1050×2450mm 79.1"×41.3"×96.5"
	Exhaust Blower	970×810×720mm 38.2"×31.2"×28.3"	1040×890×710mm 41.9"×35"×28.3"	1040×890×710mm 41.9"×35"×28.3"
Gross Wight		300Kg	340Kg	390Kg