

VOS200 Vacuum Drying Oven with Vacuum Pump

Digital Vacuum Display & Programmable Control

LAB1ST VOS200 series is an integrated vacuum drying oven engineered for precise and reliable control of both temperature and vacuum. It comes with a built-in vacuum pump and desiccant tank for simplified operation, and supports programmable vacuum and temperature cycles, making it an ideal solution for drying heat-sensitive and easily oxidized materials.

With a digital pressure sensor and microcomputer-controlled vacuum system, the VOS200 ensures efficient and repeatable drying processes.




Features

> Key Features

- **Fully Integrated Plug-and-Play Design:** Includes a built-in vacuum pump and desiccant tank for immediate operation.
- **Standard Programmable Control:** Offers comprehensive control over heating, vacuum, nitrogen purge, and time, including temperature ramp rate (slope) control for precise process management.
- **Rectangular Chamber Design:** Maximizes usable workspace volume.
- **Tempered & Bullet-Proof Double-Layer Glass Door:** Provides clear visibility into the chamber; compatible with inert gas injection.
- **Adjustable Door Seal & Integrated Silicone Gasket:** Maintains high vacuum integrity.
- **Stainless Steel Interior:** Durable, long-lasting, and easy to clean.
- **Oxygen-Free Operation:** Suitable for storage, heating, testing, and drying under inert or vacuum conditions to prevent oxidation.
- **Reduced Drying Time:** Cuts drying time by over **50%** compared to conventional vacuum drying ovens.

> Drying Process Control

- The vacuum drying process can be programmed with customizable vacuum cycles. Users can set upper and lower vacuum limits to further shorten drying times. For example:
 - Minimum Vacuum: 50Pa
 - Maximum Vacuum: 900Pa
 - Cycle Count: Up to 99 cycles
 - As cycles progress, humidity decreases rapidly, significantly accelerating drying.



> Programmable Cycle Control


- User-friendly interface with direct programming capability.
- Supports up to **15** program segments, each configurable from **1 to 9999 minutes**.
- Pre-set startup and shutdown timers for automated operation.

> Vacuum Control

- Digital vacuum control with auto-regulation via digital vacuum gauge.
- Adjustable vacuum range: **10–10⁵Pa** (minimum achievable vacuum: 133Pa).
- Control accuracy: **±1%**.
- Resistant silicon piezoresistive sensor ensures stable pressure readings, unaffected by airflow.

> Temperature Control Advantages

- **15-Segment Programmable Profiles:**
 - Standard temperature and hold-time settings ("Temperature – Dwell Time").
 - Ramp rate control using "Temperature–Time–Temperature" format for precise heating/cooling rate management.
- Full control over every stage: heating, holding, and cooling.
- **Heating Methods:** Chamber heating or shelf heating



Model	VOS200-33L	VOS200-63L	VOS200-93L	VOS200-123L	VOS200-213L	VOS200-503L	VOS200-933L
Temp. Range	RT+10°C ~ 200°C						
Display Resolution	0.1°C						
Temp. Fluctuation	±1°C						
Vacuum Range	10~10 ⁵ Pa (min. 133 Pa)						
Vacuum Sensor	Resistive Silicon Tube Pressure Sensor						
Vacuum Gauge	Digital Play						
Chamber Material	SUS304 Stainless Steel (1Cr-18Ni9Ti)						
Chamber Volume	33 L 1.17 cu ft	64 L 2.26 cu ft	91 L 3.21 cu ft	125 L 4.41 cu ft	215 L 7.59 cu ft	430 L 15.19 cu ft	910 L 32.14 cu ft
Shelves	2 pcs	3 pcs	2 pcs	3 pcs	3 pcs	4 pcs	5 pcs
Heating Method	Chamber Heating	Chamber Heating	Chamber Heating	Shelf Heating	Shelf Heating	Shelf Heating	Shelf Heating
Inner Dimensions (W×D×H, mm)	320×320×320	400×400×400	450×450×450	500×500×500	560×600×640	630×810×845	750×1160×1050
Outer Dimensions (W×D×H, mm)	550×490×1240	600×570×1390	610×590×1350	660×640×1400	720×820×1750	1000×1040×1855	1400×1395×2010
Power Input	1200 W	1800 W	1350 W	2050 W	2100 W	3800 W	5600 W
Power Supply	AC220V 50Hz	AC220V 50Hz	AC220V 50Hz	AC220V 50Hz	AC220V 50Hz	AC380V 50Hz	AC380V 50Hz

*Note: Performance data measured under no-load conditions at 20°C ambient temperature and 50% RH.

Options

Exhaust Oil Mist Filter
(recommended replacement every 6 months)

Labfirst Scientific Instruments (Shanghai) Co., Ltd.
No.248 Guanghua Road, Minhang District, Shanghai, China

CALPHA INDUSTRIES INC.
10 Corporate Park Ste 330, Irvine, California 92606, USA

LSI PROCESS EQUIPMENT CO., LTD.
701 West Georgia Street, Suite 1500, Vancouver, British Columbia, V7Y 1C6, Canada

Tel: +86-13524020331 / +1-844-452-2178

sales@lab1st.com
www.lab1st.com
eshop.lab1st.com