



LAB1ST

Semi-micro Balance



Labfirst Scientific Instruments (Shanghai) Co., Ltd.

Semi-micro Analytical Balance



Product Features

- Recipe database: 99 data can be stored
- Automatic internal calibration
- Tare/ counting/ percentage test/ animal weighing
- Baud rate is adjustable/ Density function
- Unit conversion(g/mg/ct/oz)...
- Multi-Language:English/ German/ French/ Italian/ Spanish/ Portuguese
- Keep weight values in GLP/ ISO format
- Includes date, time, and serial number
- Plastic protection cover for keypad
- Overload alarm/ fault alarm/ level indicator



Optional: Electrostatic Eliminator

Parameter

Model	PR45	PR125D	PR225D	PR225D(1)
Capacity(g)	42	62/120	82/220	102/220
Min weighing(g)		60mg		100mg
Stable Time	≤6s			
Display	LCD (white back light with black font)			
Pan size	Φ80mm			
Operation temp	5-35℃ 20-80%Humidity			
Repeat ability	0.03 / 0.1mg		0.05 / 0.1mg	0.05mg
Liner	± 0.08 / ± 0.2mg		± 0.1 / ± 0.2mg	± 0.1mg
Draft shield size	152mm×160mm×225mm			
Cal.weight	Internal Calibration			
Option	Printer			
Dimension	345×215×345mm			
G.N . weight	8.25/5.8 kg			
IONISER	Option			
Language	Italian, English, Portuguese, German, French, Spanish			
Packing Size	510x400x620mm			
Power Supply	110-230V 50/60Hz Out DC 12V 1A			

Model	PR125-ION	PR165-ION	PR215-ION
Capacity(g)	122	162	210
Min weighing(g)	100mg	100mg	100mg
Stable Time	≤6s		
Display	LCD (white back light with black font)		
Pan size	Φ80mm		
Operation temp	5-35℃ 20-80%Humidity		
Repeat ability	0.05mg	0.05mg	0.05mg
Liner	± 0.1mg	± 0.1mg	± 0.1mg
Draft shield size	152mm×160mm×225mm		
Cal.weight	Internal Calibration		
Option	Printer		
Dimension	345×215×345mm		
G.N . weight	8.3/5.85 kg		
IONISER	Standard		
Language	Italian, English, Portuguese, German, French, Spanish		
Packing Size	510x400x620mm		
Power Supply	110-230V 50/60Hz Out DC 12V 1A		