

LAB1ST

Electrochemical meter



Labfirst Scientific Instruments (Shanghai) Co., Ltd.



Model		Parameter									Channel		Security and compliance			Display Type			Data Transfer		Accessories	
Series	Model	PH	T	ORP	EC	Resistivity	TDS	ISE	DO	Ion	Single	Multi	CE	GLP	GMP	Color LCD Touchscreen	LCD Display [Color]	LCD Display [b/w]	USB	PC/RS232	Printer	Scanner
MeasureEasy	ME1	●	—	○	—	—	—	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—
MeasureEasy	ME2	●	—	○	—	—	—	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—
MeasureEasy	ME2-Pro	●	●	○	—	—	—	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—
MeasureEasy	ME3	—	—	—	●	—	—	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—
MeasureEasy	ME4	—	—	—	●	—	○	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—
MeasureEasy	ME4-Pro	—	●	—	●	—	○	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—
MeasureGo	MG10	●	●	○	—	—	—	—	—	—	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG20	●	●	○	—	—	—	—	—	—	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG30	—	●	—	●	○	○	○	—	—	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG40	—	●	—	●	○	○	○	—	—	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG50	—	●	—	—	—	—	—	●	—	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG60	○	○	○	—	—	—	—	—	○	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG70	○	○	○	—	—	—	—	—	○	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG91	●	●	○	●	○	○	○	—	○	●	—	●	●	—	—	●	—	●	●	○	—
MeasureGo	MG91-Pro	●	●	○	●	○	○	○	●	○	●	—	●	●	—	—	●	—	●	●	○	—
MeasureExcellent	M100	●	●	○	—	—	—	—	—	—	●	—	●	●	●	●	—	—	●	●	○	○
MeasureExcellent	M200	●	●	○	—	—	—	—	—	—	●	—	●	●	—	●	—	—	●	●	○	—
MeasureExcellent	M300	—	●	—	●	○	○	○	—	—	●	—	●	●	●	●	—	—	●	●	○	○
MeasureExcellent	M400	—	●	—	●	○	○	○	—	—	●	—	●	●	—	●	—	—	●	●	○	—
MeasureExcellent	M500	—	●	—	—	—	—	—	●	—	●	—	●	●	●	●	—	—	●	●	○	○
MeasureExcellent	M600	—	●	—	—	—	—	—	●	—	●	—	●	●	—	●	—	—	●	●	○	—
MeasureExcellent	M700	○	○	○	—	—	—	—	—	○	●	—	●	●	●	●	—	—	●	●	○	○
MeasureExcellent	M800	○	○	○	—	—	—	—	—	○	●	—	●	●	—	●	—	—	●	●	○	—
MeasureExcellent	M910	●	●	○	●	○	○	○	—	○	●	—	●	●	●	●	—	—	●	●	○	○
MeasureExcellent	M910-Pro	●	●	○	●	○	○	○	●	○	●	—	●	●	●	●	—	—	●	●	○	○
MeasureExcellent	M920-Multi	●	●	○	●	○	○	○	●	●	—	●	●	—	●	—	—	●	●	○	—	—

Model		Parameter										Channel		Security and compliance			Display Type			Data Transfer		Accessories	
Series	Model	PH	T	ORP	EC	Resistivity	TDS	ISE	DO	Ion	Single	Multi	CE	GLP	GMP	Color LCD Touchscreen	LCD Disply [Color]	LCD Disply [b/w]	USB	PC/RS232	Printer	Scanner	
ExplorerGo	EG10	●	—	○	—	—	—	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—	
ExplorerGo	EG20	●	●	○	—	—	—	—	—	—	●	—	—	—	—	—	—	●	—	—	—	—	
ExplorerGo	EG30	—	—	—	●	—	—	—	—	—	●	—	●	—	—	—	—	●	—	—	—	—	
ExplorerGo	EG40	—	●	—	●	—	○	—	—	—	●	—	—	—	—	—	—	●	—	—	—	—	
ExplorerGo	EG50	—	●	—	—	—	—	—	●	—	●	—	—	—	—	—	—	●	—	—	—	—	
ExplorerDirect	ED10	●	●	○	—	—	—	—	—	—	●	—	●	●	—	—	●	—	●	●	○	—	
ExplorerDirect	ED20	—	●	—	●	○	○	○	—	—	●	—	●	●	—	—	●	—	●	●	○	—	
ExplorerDirect	ED30	—	●	—	—	—	—	—	●	—	●	—	●	●	—	—	●	—	●	●	○	—	
ExplorerDirect	ED91	●	●	○	●	○	○	○	—	○	●	—	●	●	—	—	●	—	●	●	○	—	
ExplorerDirect	ED91-Pro	●	●	○	●	○	○	○	●	○	●	—	●	●	—	—	●	—	●	●	○	—	
ExplorerExcellent	E10	●	●	○	—	—	—	—	—	—	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E20	●	●	○	—	—	—	—	—	—	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E30	—	●	—	●	○	○	○	—	—	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E40	—	●	—	●	○	○	○	—	—	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E50	—	●	—	—	—	—	—	●	—	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E50-Pro	—	●	—	—	—	—	—	●	—	●	—	—	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E60	○	○	○	—	—	—	—	—	○	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E70	○	○	○	—	—	—	—	—	○	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E91	●	●	○	●	○	○	○	●	○	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E92	●	●	○	●	○	○	○	—	○	●	—	●	●	—	●	—	—	●	●	○	○	
ExplorerExcellent	E92-Pro	●	●	○	●	○	○	○	●	○	●	—	●	●	—	●	—	—	●	●	○	○	

● Standard electrode ○ Optional electrode — No such function