

**LAB1ST**



# Mixing



Labfirst Scientific Instruments (Shanghai) Co., Ltd.

# Vortex C

## Cell Disruptor

Specifications	Vortex C
Functions	Disrupts cells, yeast, bacteria, algae, fungi, etc. to prepare cell lysates for further processing
Mixing motion	High-speed circumferential vortex and oscillation function
Orbital diameter	4mm
Speed range	0-3000rpm (Adjustable)
Speed display	Scale
Capacity	8x2 mL
Power Requirements	200-240V/ 100-120V, 50/60 Hz, 60W
Dimension	127X130X160mm
Weight	2.0kg




## Features

- Durable, reliable and affordable
- Easy to operate, rapid hands-free disruption of up to 8 microtubes
- Processes at a constant, controlled shear rate
- Maintains genome integrity of the biological samples
- Facilitates subsequent purification of genomic material
- Best suited for research laboratories demanding accurate, reproducible, repeatable results


## Applications

Suitable for all sample types; Lysis of yeast and fungi; for Isolation of genomic material from faecal samples; fungi and bacteria samples; soil samples; for Isolation of PCR-Ready Genomic DNA from Soil Samples; DNA & RNA Purification; Lysis of bacteria; Lysis of soft tissues; Isolation of genomic DNA; Purification of DNA; Cell Lysis; Used to lyse cell wall of yeast; Break hard/brittle samples.


**VT1.3.1 tube adapter** Cat. No.921004  
for 48 test tubes, Ø6mm (1 ml)



**VT1.2 Tube holding rod** Cat. No.921002  
used with tube adapters




**VT1.3.2 tube adapter** Cat. No.921005  
for 15 test tubes, Ø10mm (1.5-2 mL)




**VT1.3 Universal top plate** Cat. No.921003  
Ø 100mm




**VT1.3.3 tube adapter** Cat. No.921006  
for 16 test tubes, Ø12mm (5 mL)




**VT1.3.6 Platform pad** Cat. No.921009  
Ø99mm, multiple small tubes and vessels




**VT1.3.4 tube adapter** Cat. No.921007  
for 8 test tubes, Ø16mm (10 mL)




**VT1.3.7 vacuum chuck** Cat. No.921010  
made of rubber




**VT1.3.5 tube adapter** Cat. No.921008  
for 8 test tubes, Ø20mm (10-15 mL)



**PS1.2 Microplate clamp** Cat. No.921202



**VT1.3.6 PCR tube adapter** Cat. No.921012  
for 56 test PCR tubes, strips of 8 ,  
Ø5mm (0.1mL-0.5 mL)




## Application of accessories

Accessories	Adjustable speed model (0-2500rpm)	
	Touch mode (High speed area)	Continuous mode (Low speed area)
VTI.1	Y	Y
VTI.2 + VT1.3.1	Y	
VTI.2 + VT1.3.2	Y	
VTI.2 + VT1.3.3	Y	
VTI.2 + VT1.3.4	Y	
VTI.2 + VT1.3.5	Y	
VTI.3 + VT1.3.6	Y	
VTI.3 + VT1.3.1		Y
VTI.3 + VT1.3.2		Y
VTI.3 + VT1.3.3		Y
VTI.3 + VT1.3.4		Y
VTI.3 + VT1.3.5		Y



Lab1st mixers range can be divided into different series as per mixing speed, working mode, load capacity etc:

Series	Image	Model	Speed Range	Working mode	Capacity
<b>Vortex Mixer</b>		Vortex S	0~3000rpm	Orbital (touch + continuous modes)	-
		Vortex F	3000rpm(fixed)	Orbital (touch + continuous modes)	-
		Vortex E	3000rpm (fixed)	Orbital (touch mode)	-
		Vortex C	0~3000rpm	Orbital (touch + continuous modes)	8x2mL
<b>Microplate Mixer</b>		Microster M	0~1500rpm	Orbital (continuous mode)	0.5Kg
<b>Disk Rotator</b>		Rotaster 3	10~70rpm	360° rolling (timer + continuous)	60x1.5mL, 16x15mL, 8x50mL, (available for 3 layers)
<b>Rotisserie Rotator</b>		Grillster T	10~70rpm	360° rolling (timer + continuous)	48x1.5mL, 24x15mL, 24x50mL
<b>Disk Rotator</b>		Rotaster 1 basic	0~80rpm	360° rolling (continuous)	48x2mL, 16x15mL, 8x50mL
<b>Rotisserie Rotator</b>		Grillster basic	0~80rpm	360° rolling (continuous)	32x2mL, 20x15mL, 16x50mL
<b>Rocker</b>		See-saw 30 basic	0~30rpm	30±3° rocking (continuous)	10x10mL, 4x50mL
		See-saw 30L	0~30rpm	30±3° rocking (continuous)	16x10mL 6x50mL
		See-saw 30D	0~30rpm	40±3° rocking (continuous)	32x10mL, 12x50mL
<b>Roller</b>		Rollster 6	0~70rpm	360° rocking&rolling (timer + continuous)	4Kg
		Rollster 6 basic	0~70rpm	360° rocking&rolling (continuous)	4Kg

Lab1st shaker series can be divided to different types as per different motion models:

Series	Image	Model	Plate Size	Display	Speed Range	Max. Load Capacity
<b>Orbital Shaker</b>		Orbiter 330	32x32cm	LCD	100-500rpm	7.5Kg
		Orbiter 250	25x24.5cm	LCD	100-800rpm	2.5Kg
		Orbiter 330 basic	33.5x33.5cm	LCD	70-400rpm	3Kg
		Orbiter 250 basic	26.8X26.8cm	LED	40-200rpm	2Kg
<b>Linear Shaker</b>		Recipster 330	32x32cm	LCD	100-350rpm	7.5Kg
		Recipster 250	25x24.5cm	LCD	100-350rpm	2.5Kg
<b>Rocking Shaker</b>		Wavester 330	30x28.5cm	LCD	10-70rpm	10Kg
		Wavester 250 basic	26.8x26.8cm	LED	10-80rpm	2Kg
<b>3D Shaker</b>		3D-Ster 330	30x28.5cm	LCD	10-70rpm	5Kg